wellion[®] SENSOR

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Applying a sensor 70

Removing a sensor 73

Wellion AiDEX SENSOR APP – Diabetes Management System 74

Caring for your Continuous Glucose Monitoring System (CGMS) 74

Specifications **75**

Electromagnetic compatibility **77**

Symbols **81**

Important safety information

Indications for use

The Wellion SENSOR Continuous Glucose Monitoring System (CGMS) is indicated for continuous or regular monitoring of glucose levels in subcutaneous tissue, and is used for daily detection and self-management of glucose levels in people 14 years or older. It is intended for use by users at home and in healthcare facilities.

Interpretation of the Wellion SENSOR results should be based on the glucose trends and several sequential readings over time. The Wellion SENSOR CGMS also aids in the detection of episodes of hyperglycemia and hypoglycemia.

For any changes in therapy, do not only use the glucose values of the Wellion SENSOR CGMS, but also the comparative values of a blood glucose meter.

Users

- The Continuous Glucose Monitoring System is suitable for users with diabetes requiring regular blood glucose monitoring.

Contraindications

Users who:

- Suffer from alcoholism, drug abuse, severe mental disorders (such

as depression, schizophrenia).

- Are unconscious.
- Are unable to understand or master the usage of the device.
- Have severe hearing or vision impairment.
- User below the age of 14 or are not capable of administering diabetes therapy themselves or who always need to use Continuous Glucose Monitor under strict supervision of parents or qualified healthcare professionals.

Precautions

- CGMS readings should only be used as a reference for the supplemental monitoring of diabetes mellitus and should not be used as a basis for clinical diagnosis.
- The CGMS should be completely removed before Magnetic Resonance Imaging (MRI).
- The CGMS contains many small parts that can be dangerous if swallowed.
- During rapid changes in blood glucose (more than 0.1mmol/L per minute), glucose levels measured in interstitial fluid by the CGMS may not be the same as blood glucose levels. When blood glucose levels drop rapidly, the sensor may produce a higher reading than the blood glucose level; Conversely, when blood glucose levels rise rapidly, the sensor may produce a lower reading than the blood glucose level. When experiencing symptoms of hypoglycemic or hyperglycemic, the sensor readings are checked by using a blood

glucose meter for a fingertip blood glucose test.

- Severe dehydration or excessive loss of water may result in inaccurate results. When you suspect you are dehydrated, consult a healthcare professional immediately.
- If you think the CGMS sensor reading is inaccurate or inconsistent, use a blood glucose meter to test your blood glucose level or calibrate the glucose sensor. If the problem persists, remove and replace the sensor.
- While extensive user testing was done on Wellion SENSOR CGMS in Type I and Type II diabetics, the study groups did not include women with gestational diabetes.
- The performance of the CGMS has not been evaluated when used with another implantable medical device, such as a pacemaker.
- Only Wellion SENSOR consumables should be used with the CGMS.
- If the product is not working properly or has been damaged, stop using the product.
- When used in a medical facility, the operator should wear gloves while using.

Product components


Name: Wellion SENSOR continuous glucose monitoring system

The glucose sensor system consists of the following components:

Applicator with Built-in Sensor

Transmitter

Receiver

Portable Diabetes Assistant (PDA)

Wellion AiDEX
SENSOR App*

*You can download the Wellion AiDEX SENSOR App from App Store (iOS) or Google Play Store (Android)

Component description


Wellion SENSOR Portable Diabetes Assistant (PDA)

PDA may be purchased separately – for more information about usage refer to the enclosed manual

EN

Wellion SENSOR Transmitter


The transmitter is used together with the glucose sensor base. When the sensor and transmitter are activated on one's body, it measures and stores glucose readings as well as transmitting the readings to the PDA/the Wellion AiDEX SENSOR mobile App.

Wellion SENSOR CGMS Pack (=sensor base)

This package contains a sterilized sensor applicator with integrated sensor base.

The package is sterilized by gamma irradiation.

Using a new glucose sensor

Applying a sensor

- Recommended areas for sensor application includes the abdomen, the outside and the back of the upper arm. For best performance, avoid excessive motion which may weaken the sensor and its adhesive tape. Avoid accidental knocking off the sensor. It is strongly recommended to put sensor on your abdomen rather than arms or other areas.


- Before application, use alcohol pads to disinfect your selected skin area and wait for a minute to let the skin dry. Change your sensor sites regularly to avoid discomfort or skin irritation.
- Open the sensor package.

Note: Ensure the sensor is not expired by checking its expiration date. If the sensor is expired or the sterilized package is damaged, please do not use it


- Take out the sensor applicator by turning it anti-clock wise.
- Place the applicator on the top of the desired sensor site. Press it firmly against skin and press the button to launch the sensor. After inserting the sensor, wait a few seconds to allow the sensor base patch to adhere to the skin.
- Remove the applicator. Now, the sensor should be successfully attached.

EN


Note: Applying a sensor may cause bruising or bleeding. If bleeding does not stop, remove the sensor and apply a new sensor to a different location.

Release button

- Align the transmitter to the sensor base and press the transmitter firmly onto the sensor base until it snaps onto the sensor base and attaches to it securely.

Note: If the orientation of the transmitter and sensor base mismatch, the transmitter cannot be installed properly.


- The recyclable applicator is meant for single use only. Dispose the used sensor package and sensor applicator according to your local regulations.
- Follow the steps described above carefully and be sure to only use components that are made by MicroTech Medical. Use of unapproved components can result in injury.

EN

Removing a sensor

- Remove the old sensor when the PDA/Wellion AiDEX SENSOR app indicates that the sensor has expired or if you feel any site irritation or discomfort.
- Carefully pull the tape that holds the sensor to the skin from the edge and slowly uncover it until the entire sensor is removed. You can use warm soapy water to remove any remaining sticky residue.
- Press the two locking arms on one side of the sensor base, pull the transmitter away from the sensor base, and then remove the transmitter so that it can be used again.
- The sensor is meant for one-time use only. Dispose the used sensor according to your local regulations.

Wellion AiDEX SENSOR APP – Diabetes Management System

The Wellion AiDEX SENSOR App is used combined with the Wellion SENSOR CGMS for continuous glucose monitoring.

EN For more information visit www.wellion.at

Caring for your Continuous Glucose Monitoring System (CGMS)

Cleaning the transmitter

- Disinfect the transmitter with an alcohol wipe.
- Do not use solvents, nail polish remover, or paint thinner to scrub the outer surface.
- Keep the transmitter dry.
- Do not use any lubricant.

Note: Do not immerse the transmitter in water or other liquids.

Disposal

Dispose transmitters and sensors that you do not use any more in compliance with local regulations for electronic devices, batteries, sharps, and biohazard materials. Please do not discard old products or accessories directly. For further information on how to dispose of system components, please contact customer service.

Transportation

Avoid placing heavy weight on top of the PDA and transmitter. Avoid direct sunlight and rain.

Storage

If you are temporarily not using the transmitter, or sensor system, store the components in a cool, dry, clean and well-ventilated area.

The transmitter and sensor are precision instruments. If they fail, they can only be returned to the manufacturer to repair. No third-party individuals or organizations are allowed to perform repairs. Circuit diagrams and component lists are not provided in the manual.

Specifications

Item		Subcomponent		
		Transmitter	Sensor	PDA
14 days	Model number & Use Life	WELL19-S14D00T (G7-T01)	WELL19-S14D (G7-S01)	WELL19-S00PDA WELL19-S00PDAMM (G7-P02)
		4 years	14 days	4 years
Durability: 1 year				
10 days	Model number & Use Life	WELL19-S10D00T (G7-T01A)	WELL19-S10D (G7-S01A)	WELL19-S00PDA WELL19-S00PDAMM (G7-P02)
		4 years	10 days	4 years
Durability: 1 year				

Operating temperature	5-40°C (41-104°F)		
Operating humidity	10-93% (non-condensing)		
Storage and transportation temperature	-20 - 60°C	4 - 30°C	-20 - 60°C
Storage and transportation humidity	5-95% (non-condensing)		
Storage and transportation pressure	700hpa-1060hpa		
Ingress protection level	IPX7		IPX0
Detection range	2,0 mmol/L-25,0 mmol/L (36 mg/dl - 450 mg/dl)		
Wireless frequency and bandwidth	Frequency: 2,402GHz-2,48 GHz Bandwidth: 1Mbps		
Wireless modulation	GFSK		
Radiated power	-2dBm		

Electromagnetic compatibility

These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.

EN

Portable and mobile RF communication interference may have an impact on the device.

The use of accessories other than those specified for the device is not recommended. They may result in increased emissions or decreased immunity of the device.

The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

The basic performance is described in the following table:

Performance	Description
Measurement accuracy	When the glucose concentration is >4.2 mmol/L (75mg/dL), the sensor accuracy deviation does not exceed $\pm 20\%$; When the glucose concentration is ≤ 4.2 mmol/L (75mg/dL), the accuracy deviation does not exceed $\pm 1,1$ mmol/l (20mg/dL)

IEC60601-1-2:

Emissions test	Compliance
Radio Frequency emissions CISPR 11	Group 1, Class B

IEC 60601-1-2:

Guidance and manufacturer's declaration – electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge(ESD) IEC 60601-4-2	± 8 kV Contact ± 15 kV Air	± 8 kV Contact ± 15 kV Air
Electrical fast transient burst IEC 61000-4-4	± 2 kV power cord ± 1 kV input/output	± 2 kV power cord ± 1 kV input/output
Surge IEC 61000-4-5	± 0.5 kV and ± 1kV line(s) to line(s).	± 0.5 kV and ± 1 kV line(s) to line(s).
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4- 11	0% 230 V for 1 cycle 0% 230 V for 0.5 cycle at 8 phase angles 70% 230 V (30% dip in 230V) for 25/30 cycles 0% 230V for 250/300 cycles	0% 230 V for 1 cycle 0% 230 V for 0.5 cycle at 8 phase angles 70% 230 V (30% dip in 230V) for 25/30 cycles 0% 230V for 250/300 cycles
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m

IEC 60601-1-2:

Guidance and manufacturer's declaration – electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-67	6 Vrms 150 kHz - 80 MHz	6 Vrms 150 kHz - 80 MHz
Radiated RF IEC 61000-4-3	10 V/m 80MHz-2,7GHz	10 V/m 80MHz-2,7GHz

Recommended separation distances between portable and mobile RF communications equipment and the device

EN These devices are intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment

Maximum rated output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150kHz-80MHz $d=1.2 \sqrt{P}$	80 MHz-800 MHz $d=1.2 \sqrt{P}$	800 MHz-2.5 GHz $d=2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Symbols

Class II equipment

Single use only

Type BF applied part

Storage temperature

Sterilized by radiation

Refer to instruction manual/booklet

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MicroTech Medical (Hangzhou) Co., Ltd.
LOTUS NL B.V.