w wellion®

LIPID *pro*

TESTING SYSTEM

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CONTENT

Introduction	40
Getting Started	41
Components	43
Meter	43
Test Strips	46
Control Devices	47
Initial Setup	48
Turn on the meter	48
Coding the meter	49
Meter Setup and Options	49
Test Number Setup	50
Optical Check	50
Specimen Type Set Up	51
System Setup	52
Date Setup	52
Time Setup	53
Test Number Reset Setup	53
Unit Setup	54
Sound Setup	54
Delete Data Setup	54
Testing	55
Specimen Collection	55
Test Processing	57
Data/Communication	59
Deleting Data	60
Memory/Database	60

Optical System Check	61
Quality Control	62
Maintenance	62
General Cleaning	62
Disinfection Process	63
Replacing the Batteries	64
Precautions	65
Troubleshooting	66
Meter Specifications	68
Index of symbols	69
Warranty	70

The Wellion LIPID pro TESTING SYSTEM is intended for the quantitative determination of Total Cholesterol (TC), High Density Lipoprotein Cholesterol (HDL), Triglycerides (TG), the calculated ratio of TC/HDL and Low Density Lipoprotein Cholesterol (LDL) in capillary and venous human whole blood, plasma and serum. The system consists of a portable meter that analyzes the intensity and color of the reagent area of a test strip, ensuring quick and accurate results.

The Wellion LIPID pro TESTING SYSTEM provides results in less than 2 minutes. The meter can store up to 200 results and records can be transferred to a computer for further analysis using the USB port. The meter can be operated by 3 AAA (1.5V) batteries.

To ensure accurate results:

- Read instructions carefully and complete any necessary training before use.
- Use the code chip that is included in each box of test strips.
- Only use the Wellion LIPID 5in1 TEST KITS with the Wellion LIPID pro TESTING SYSTEM.
- For in vitro diagnostic use only.
- For professional use only.
- Sample type: Fresh capillary blood, heparinized or EDTA venous whole blood, serum and heparinized plasma can be tested.
- Please perform tests within the expiry date of all components.

Before testing, read the instructions carefully and learn about all the components of the Wellion LIPID pro TESTING SYSTEM. Depending on the package type, some of the components may need to be purchased separately. Please check the list of contents on the outer box for details on which components are included with your purchase. The following items are needed to perform a test:

Included in the Wellion LIPID pro TESTING SYSTEM set:

Wellion LIPID Pro
TESTING SYSTEM

3x AAA batteries

Control Strip

Separately available:
Wellion SafetyLancets

Included in the Wellion LIPID 5in1 TEST KIT (separately available):

Capillary Transfer Tube/
Dropper

Front Back

Test Strip


Code Chip

- **Wellion LIPID pro TESTING SYSTEM:** Reads the test strips and displays the concentrations of TC, HDL, TG, and calculated LDL and TC/HDL values.
- **Test Devices:** Part of the TEST KIT, work with the meter to measure the concentrations of TC, HDL, TG and calculated LDL and TC/HDL values.
- **Code Chip:** Automatically calibrates the meter when inserted into the meter.
- **Capillary Transfer Tubes/Droppers:** Collects capillary blood from fingertip for accurate results (25 μ L).
- **AAA Batteries:** Provides mobility for testing.
- **Carrying Case:** Allows you to transport the device.
- **User's Manual:** Provides detailed instructions on using the Wellion LIPID pro TESTING SYSTEM
- **Quick Reference Guide:** Provides a brief overview of the Wellion LIPID pro TESTING SYSTEM and its testing procedures.
- **Test Strip Package Insert:** Provides detailed instructions on using the Wellion LIPID 5in1 TEST KIT.
- **Safety Lancets:** Used to draw blood specimens. Discard after use.
- **Control Strip:** Verifies the proper operation of the meter by checking that the meter can detect a pre-calibrated value.


The Wellion LIPID pro TESTING SYSTEM reads the test strips and displays the concentrations of TC, HDL, TG, the calculated value of LDL and the ratio of TC/HDL. Use the below diagram to become familiar with all the parts of the meter.

Meter

* the printer is sold separately

During testing, the Wellion LIPID pro TESTING SYSTEM will display icons showing the status, options available, and prompts for testing:

- | | |
|-----------------------|-----------------------|
| 1. Code | 8. Battery |
| 2. SET option | 9. Sample type |
| 3. Test item | 10. Measurement units |
| 4. Memory | 11. Test result area |
| 5. Blood drop symbol | 12. Sound icon |
| 6. Date | 13. Test number |
| 7. Test device symbol | |

Code: Shows the code number of the test strips.

SET option: Setup the meter functions.

Test Item: Shows which item is being tested.

Memory: Indicates a test result is being recalled from memory.

Test Strip and Blood Drop Symbols: Indicates when to insert test strip or apply specimen.

Date: Shows the current date or date tested.

Battery: Appears when the battery should be replaced.

Sample type: Displays the type of blood.

Measurement Units: Displays the units of the test result.

Test Result Area: Displays the test result or menu options.

Sound Icon: Appears when the sound is turned on.

Test Number: Indicates the assigned test number.

- Avoid water or other liquids on or inside the meter.
- Keep the Strip Channel clean.
- Keep the meter dry and avoid exposing it to extreme temperatures and humidity.
- Do not drop the meter or get it wet. If the meter is dropped or has gotten wet, ensure the meter is working properly by running an optical check. (Refer to Optical System Check for details)
- Do not take the meter apart. Taking the meter apart will void the warranty.
- Refer to Maintenance for details on cleaning the meter.
- Keep the meter and all associated parts out of reach of children.

Note

Follow proper precautions and all local regulations when disposing of the meter and used batteries.

All Lipid Testing Systems Preventive Warnings with Regard to Electromagnetic Compatibility (EMC)

1. This instrument is tested for immunity to electrostatic discharge as specified in IEC 61000-4-2. However, use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets, etc.) may cause erroneous results.
2. These instrument and immunity requirements are described in EN 61326-1 and EN 61326-2-6. Do not use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with proper operation of the meter.
3. For professional use, the electromagnetic environment should be evaluated prior to operation of this device.

Test Strips

The test strips of the Wellion LIPID 5in1 TEST KIT work with the Wellion LIPID pro TESTING SYSTEM to measure the lipid concentration in whole blood, plasma and serum.

The test strips of the Wellion LIPID 5in1 TEST KIT can detect TC, HDL and TG values at the same time. The ratio of TC/HDL and the value of LDL are calculated by the meter additionally.

Insert Arrow: Located on the front of the test strip, the arrows indicate the direction in which the test strip should be inserted into the meter.

Specimen Application Area: After the strip is inserted into the Strip Channel, apply the correct specimen volume (25 μL) to the region in the center of the test strip.

Handle: Located at the end of the test strip, the handle is used to insert and remove the test strip from the meter.

Test Area: Located at the back of the test strip. The meter will detect and read this area to give results of lipid levels.

Specimen Application

For best results, fill the Specimen Application Area with the correct specimen volume (25 μL). Incorrect results may occur, if the specimen is not applied correctly or if the Specimen Application Area is not filled with the correct amount.

After applying the specimen, ensure that the Specimen Application Area is completely covered. The Specimen Application Area should remain covered throughout the entire test. If the Specimen Application Area is not covered or if there is too much specimen covering the Specimen Application Area, repeat the test with a new test strip.

Note

If the specimen applied to the Specimen Application Area is not enough, do not add more specimen to the test strip. Instead, retest with a new strip. If the E-5 Error or another error appears on the display, please discard the used test strip and retest with a new strip.

Code number

Printed on each package of test strips is a code number **CODE**, lot number **LOT**, unopened expiration date , and test strip quantity .

Test Strip Precautions and Instructions for Use

- The test strips of the Wellion LIPID 5in1 TEST KIT should be stored in their canister to keep them in working condition.
- Do not store test strips outside of their package. Test strips must be stored in the original package.
- Do not transfer test strips to a new package or any other container.
- For in vitro diagnostic use. Test strips are to be used only outside the body for testing purposes.
- Use the test strip immediately after removing it from the vial.
- Do not use test strips that are torn, bent, or damaged in anyway. Do not reuse test strips.
- Before performing a test, make sure that the code number on the meter display matches the number shown on the test strip vial and on the ink-jet printing on the code chip.
- Refer to the test strip package insert for more details.

Control Devices

The Wellion LIPID Control Strips are devices containing a reference pad, which work with the Wellion LIPID pro TESTING SYSTEM to ensure the optical system is working properly. After the control strip is inserted into the meter, the meter's optical system detects the color intensity of the control strip. The meter displays YES or no to indicate whether the meter is functioning properly. Refer to Optical System Check for details.

Precautions

- Store in the closed canister at a temperature within 2-30°C (36-86°F). Avoid exposure to direct sunlight, extreme temperatures, and humidity.
- Control strips should be stored in their tightly capped vial to keep them in working condition.
- Do not freeze.
- Keep the control strips clean. Do not touch the test area of the strip.
- Remove the control strip for immediate out of the vial. Put the control strip back and close the vial immediately after use. Do not use contaminated, discolored, or damaged control strips.
- Do not use after the expiration date.
- For in vitro diagnostic use only.

Storage and Handling

- Store test strips in a cool, dry place. Store away from heat and direct sunlight.
- Transport and store in its closed canister within 2-30°C (36-86°F) temperature and less than 90% humidity.
- Do not freeze.
- Replace the cap on the strips vial immediately after removing a strip. Expired strips may cause incorrect test results.

Note

The expiration date is printed in a Year-Month format. For example, 2024-01 is January, 2024.

INITIAL SETUP

Before testing, ensure the following procedures are followed.

Turn on the meter

The meter can be operated using 3 AAA batteries (1.5V).

To use the meter with batteries, insert 3 AAA batteries (1.5V) into the battery compartment on the back of the meter.




Press to turn the meter on after the batteries are inserted. The screen will briefly display all of the LCD symbols. Observe the LCD at startup to ensure all segments and display elements are turned on. After the power-on diagnostic check, the Initial Screen will be displayed. The meter will automatically turn off after 5 minutes of inactivity.

Coding the meter

Each time a new box of test strips is used, the new code chip included in the box must be inserted into the meter. Compare the code number on the code chip from the box with the code number printed on the test strip canister label. Results may be inaccurate if the two numbers are not identical. Insert the new code chip into the code chip slot of the meter.

The code chip should remain in the meter. Do not take it out until a new box of test strips is needed. The code number will appear on the Initial Screen after startup.

If the code chip is not properly inserted into the code chip slot or if it is missing, the meter will display three dashes as shown below

METER SETUP AND OPTIONS

With the meter turned on, press SET button to enter the Meter Setup mode.

Press ◀ or ▶ to display several setup sub-modes:

n.SET	Test number setup. The test number can be set from 1 to 999.
--------------	--


CHE	Optical Check mode. Refer to Optical System Check.
SEt SE BL	Specimen type setup.
SYS	System setup, including date, time, test number reset, units, and sound.
dEL	Memory Delete mode. Refer to Data/Communication.
EIt	Exit setup modes and save changes when ϕ is pressed. The meter will automatically return to the Initial Screen.

Press ϕ to enter the mode when the desired sub-mode is displayed.

Test Number Setup

From the **n.SET** screen, press ϕ to enter Test Number Setup. The test number can be set to any number from 1 to 999.

Press \blacktriangleleft or \blacktriangleright until the correct test number is displayed. To quickly cycle to the desired test number press and hold one of the arrow button arrow buttons. Press ϕ to save and return to the Meter Setup screen.

Optical Check

From the **n.SET** screen, press \blacktriangleright to enter the **CHE** Setup.

Note

The control strip is intended for checking the optical system. Press ϕ to enter Check control strip mode. Insert a control strip into the strip channel in the same direction as the arrows indicate on the strip. Ensure that the test strip is inserted all the way, then press ϕ to confirm.

If the meter displays YES, the meter is normal. If the meter displays no, the meter is abnormal, as shown below.


If the meter displays no, check if the control strip is contaminated, bent or damaged. If there are any visible signs of damage or contamination, discard the control strip and retest a new control strip.

Press the **⏪** button to return to the Setup Screen.

Specimen Type Set Up

From the **n.SET** screen, press **▶** to enter the **SET SE BL** screen. Press **⏪** to set specimen type.

Press **◀** or **▶** to set specimen type to either **BL** or **SE**. When specimen type is set to **BL**, control solution, fresh capillary blood, EDTA or heparinized venous whole blood can be used. When specimen type is set to **SE**, serum and heparinized plasma can be used. Press **⏪** to save and return to the setup screen.


Note

BL means whole blood, **SE** means plasma and serum

System Setup

From the **n.SET** screen, press **▶** until **SYS** Setup.

Press **⏻** to start system setup. The first option sets the clock to either 12 or 24 hour mode. Press **◀** or **▶** to switch between the two settings.

Date Setup

Press **⏻** to save and move to Date Setup Y-M-D (Year-Month-Day).




The year will flash at the bottom of the display with Y indicating year setup. Press **◀** or **▶** until the correct year is displayed.

Press **⏻** to save and move to the Month and Day Setup.

The date will appear at the bottom of the display separated by a single dash (-), with the month flashing. M will also appear indicating month setup. Press **◀** or **▶** until the correct month is displayed.

Press **⏻** to save. The day will flash and **d** will appear indicating day setup. Press **◀** or **▶** until the correct day is displayed.


Press **⏻** to save and proceed to Time Setup.

Time Setup

Hours and minutes will appear at the bottom of the display separated by a colon, with the hour flashing.


Press ◀ or ▶ until the correct hour is displayed.
Press ⏻ to save and move to Minutes

Note


The meter will display AM or PM if the 12H time setting is chosen.

Minutes will flash. Press ◀ or ▶ until the correct Minutes are displayed. Press ⏻ to save and move to Test Number Reset Setup.


Test Number Reset Setup

Press ◀ or ▶ to turn the test number reset On or OFF. The test number will reset to 1 for each new day of testing when the test number reset is turned on.


or

Press ⏻ to save and move to Unit Setup.

Unit Setup


Set the units to either mg/dL or mmol/L. Press ◀ or ▶ to switch between the two settings.

Press ⏻ to save and move to Sound Setup.

Sound Setup

Press ◀ or ▶ to switch the sound either On or OFF. The Sound Symbol will appear on the display when the sound is turned on.


Press ⏻ to save and return to **n.SET**.

Delete Data Setup

Press ▶ to enter the dEL Setup.

Press ⏻ to select dEL Setup.

Press ⏻ again to delete the MEM data.

It will automatically return to **n.SET** when the MEM data is deleted.

Before performing any test, the user should review the Wellion LIPID pro TESTING SYSTEM Manual for detailed instructions. The following steps show how to use each component to measure the lipid concentration.

Specimen Collection

1. Use fresh capillary blood from the fingertip.
2. Use heparinized or EDTA venous whole blood, serum and heparinized plasma specimens. Please refer to Professional Testing below.

Note

Before testing, choose a clean, dry work surface. Review the procedure and make sure all of the items needed to obtain a sufficient amount of blood are available.

Testing with heparinized or EDTA venous whole blood, serum and heparinized plasma

For heparinized or EDTA venous whole blood, serum and heparinized plasma, mix the specimen well, then collect specimen (25 μ L) into a plastic/glass capillary transfer tubes or pipette. Apply it to the center region of the Specimen Application Area of the device. Do not touch the test devices with the pipette or tube.

- Specimen must be tested within 8 hours of collection.
- Mix the specimens well before testing, in order to ensure the cellular components are evenly distributed.
- Allow the specimen to come to operating temperature (15-40°C or 59-104°F) for approximately 15 minutes, if the specimen has been refrigerated.
- Anticoagulants other than EDTA and heparin are not recommended.


Note

Refer to NCCLS Documents H3-A6, Collection of Diagnostic Blood Specimens by Venipuncture.

Testing with fingertip blood

Wipe away the first drop of blood. Apply light pressure to obtain a second drop of blood. Collect capillary blood (25 μ L) using a Capillary Transfer Tube or pipette.

For use with the Capillary Transfer Tube, hold the tube slightly downward and touch the tip of the Capillary Transfer Tube to the blood specimen. Capillary action will automatically draw the specimen to the fill line and stop.

Note

The Capillary Transfer Tube will fill automatically. Make sure that blood covers the black line, or it will be difficult to squeeze the blood out. Never squeeze the Capillary Transfer Tube while sampling.

Align the tip of the Capillary Transfer Tube with the center hole of the Specimen Application Area of the test strips to apply the second drop of blood (25 μ L).

Note


Do not touch the test strip with the Capillary Transfer Tube or pipette. The capillary blood should be tested immediately after collection. The use of a Capillary Transfer Tube or pipette is recommended for accurate results.

Blood specimens can be obtained by using a Wellion SafetyLancet.


Collect the blood sample

Prior to testing, make sure the patient's hand is warm and relaxed before collecting the capillary blood specimen. Use warm water to increase blood flow if necessary. Massage the hand from the wrist up to the fingertip a few times to encourage blood flow.

Clean the testing site with an alcohol swab or by washing the hands with warm soapy water and then dry the testing site thoroughly.

Put the safety lancet against the fingertip to prick the finger. Gently massage from the base to the tip of the finger to obtain the required blood volume. Avoid smearing the drop of blood. Rotation of sites is recommended. Repeated punctures in the same spot can make the fingers sore and callused.

Note


Don't use an infection swab containing iodine. This can give inaccurate results.

Note


For professional use, please refer to NCCLS Documents H04-A6, Collection of Diagnostic Capillary Blood Specimens.

Test Processing

Ensure the meter is set up properly, as described in previous sections. Turn the meter on. The screen will briefly display all the LCD symbols. Observe the LCD at startup to ensure all segments and display elements are turned on. There should be no missing icons or elements. The meter will briefly show a blank display.


After startup, the Initial Screen will be displayed. Ensure the code chip is inserted. Compare the number showed on the display with the code number printed on the foil pouch. Refer to Initial Setup. The test strip symbol will flash when the meter is ready for the strip to be inserted.


Check that the specimen type displayed on the meter LCD is the same as the specimen type tested. If not, set the correct specimen type. Refer to Specimen Type Set Up.

Testing

Insert a test strip into the Device Channel in the same direction as the arrows indicate on the test strip.


Ensure that the test strip is inserted to the end of the Strip Channel. The blood drop symbol will flash when the meter is ready for the specimen to be applied.


Apply the blood specimen (25 μL) to the center of the Specimen Application Area on the test strip.

Note

For testing capillary blood, use the second drop of blood for accurate results. The meter will begin testing automatically with three dashes in a line flashing on the display indicating that the test is in progress.

Results will be displayed within 2 minutes. Press ► to view the results.

Remove the used test strip. The meter will return to the Initial Screen and is ready to perform a new test.

Note

Discard all blood specimens, used test strips, and materials carefully.

Follow proper precautions and obey all local regulations when discarding blood specimens and materials.

Perform daily cleaning when testing is completed for the day.

Refer to the Maintenance section.

The meter will automatically turn off after 5 minutes of inactivity or when ϕ is pressed. Remove the batteries if the meter will not be used for an extended period.

DATA/COMMUNICATION

Note

Up to 200 test records are automatically stored in the memory. After 200 test records are stored, the oldest test record will be replaced by a new record. For example, if 200 records are stored in the memory, the next test result (201) will replace the first result in the memory.

Deleting Data

To delete all data from the meter database, enter the Setup Menu. Refer to Meter Setup and Options for more details. Press ◀ or ▶ until dEL is displayed.

Press ⏻ to enable data deletion, MEM will be displayed.

Press ⏻ until the meter returns to the Setup Menu.




Memory/Database

From the initial test screen, press ◀ or ▶ to proceed to the memory/database.

Press ◀ or ▶ to view the memory of the tests. For all Wellion 5-in-1 test strip results, select the option TC HDL TG. Press ⏻ to enter the selected memory screen. The screen will show the latest results. Press ◀ or ▶ to choose the test number and view each record in the date/time sequence.

To view the test results, press ⏻ to proceed to the record. Then press ◀ or ▶ to view results of TC, HDL, TG, TC/HDL and LDL. Press **SET** button to exit.

If no data is stored, the meter will display one dash (-) and MEM. Press SET button or ⏻ to return to the testing screen.

OPTICAL SYSTEM CHECK

Press ◀ or ▶ from the Setup Screen to select the Optical Check mode, as shown below.

Note

- The control strip is intended to check the optical system.
- Allow the control strips and the meter to reach operating temperature (15-40°C or 59-104°F) prior to testing.
- The optical check should be performed under normal lab lighting conditions. Do not perform under sunlight or extreme lighting conditions.

Press ⏻ to enter this mode. The meter will flash the test strip symbol as shown on the picture.

Insert a control strip into the Strip Channel. Follow the direction of the arrows indicated on the strip. Ensure that the control strip is inserted all the way.

Press ⏻ to start the optical check. If the meter displays YES, the meter is normal. If the meter displays no, the meter is not functioning properly.

or

If the meter displays no, check the control strip for contamination or damage. If there are any visible signs of damage or contamination, discard the control strip and retest using a new strip. Press ⏻ to return to the Setup Screen.

QUALITY CONTROL

Each lab should use its own standards and procedures for performance. Test known specimens/controls at each of the following events in accordance with local, state, and/or federal regulations or accreditation requirements:

- Each new day of testing
- When a new package of test strips is opened
- When a new operator uses the meter
- When test results seem inaccurate
- After performing maintenance or service on the meter

If QC tests do not provide expected results, perform the following checks:

- Ensure that the test strips used are not expired.
- Ensure that the test strips are fresh from a new package.
- Ensure that the control strips are not expired.
- Repeat the test to ensure no errors were made during the test.

MAINTENANCE

Proper maintenance is recommended for best results.

General Cleaning

For best results, the meter should be cleaned after each day of testing.

Meter Surface

A cotton cloth can be used to clean the surface of the meter. Use a damp cotton cloth if necessary.

A dry, soft cloth may be used to clean the LCD and the sensor area. It is recommended that the meter is stored in the carrying case after each use. Avoid getting liquids, residue, or control solutions in the meter through the Strip Channel, Code Chip Slot, or USB Port.

Test Strip Holder

Remove the Test Strip Holder by pressing in on the middle of the Test Strip Holder and sliding it out from the meter. Wipe it down with a damp cloth or a mild detergent. Dry it with a dry, soft cloth. Slide the Test Strip Holder back into the meter by laying it flat on the meter. Firmly press down on the two sides of the Test Strip Holder with your thumb and push it in until it clicks into place.


Note

Do not use organic solvents, such as gasoline or paint thinner. This will cause damage to the meter.

Meter Sensor Area

Remove the Test Strip Holder as described in the previous section. Wipe down the Meter Sensor Area with a cotton swab. Do not scratch the transparent window covering the sensors.

Note

Do not use bleach or alcohol to clean the Meter Sensor Area. This will cause damage to the meter.

Disinfection Process

The disinfection process should be performed before each test to prevent potential infectious disease Transmissions through blood-borne pathogens.

Cleaning Before Disinfection and How to Disinfect

Before disinfection, use EPA Registered towelette/wipes with active ingredients of Isopropyl alcohol to clean the meter. Use these towelette/wipes to remove any stains/debris. The cleaning before disinfection ensures stains or debris are removed before disinfection for an effective sterilization.

For disinfection, please use a fresh EPA Registered towelette/wipe with active Isopropyl alcohol to wipe the meter. Be sure to wet the entire outer meter surface thoroughly. The outer meter surface must remain visibly wet for one full minute. After wiping, allow the meter to air dry completely before using the meter again.

Note

Avoid inserting the towelette/wipe into the inside of the Code Chip Slot and the USB Port when performing cleaning before and during disinfection.

Disinfection Frequency

The meter disinfection process should be performed for the time of use. This ensures that your meter will operate properly with regular disinfection during the meter's life. Check normal meter electronic operations regularly. Do this by ensuring the LCD display shows all segments once the meter is turned on before testing.

Replacing the Batteries

When the battery icon  is flashing, the batteries are low and should be replaced as soon as possible. An E-4 error message will appear if the batteries are too low to perform any more tests. The meter will not function until the batteries are replaced.

Make sure the meter is off before removing the batteries. Turn the meter over to locate the battery cover. Press the battery cover tab on the top and lift the cover to open it. Remove and discard the old batteries. Insert three new AAA batteries into the battery compartment, alternating orientation up and down as indicated in the bottom of the battery compartment.


Close the battery cover and make sure that it snaps shut. Recheck and reset the clock setting if necessary, after replacing the batteries to ensure that the time is correctly set. Refer to Initial Setup.

Note

Do not discard batteries with household waste. Follow local regulations for disposal.

Follow the precautions listed below to ensure accurate results and proper operation of the meter.

- The protection provided by the equipment may be impaired if used in a manner not defined in this instruction manual.
- Wear gloves to avoid contact with potentially hazardous biological specimens during testing.
- Avoid storing or operating the meter in direct sunlight, excessive temperatures, or high humidity. Refer to Appendix 1 Meter Specifications for operating condition requirements.
- Keep the unit clean. Wipe it frequently with a soft, clean, and dry cloth. Use a damp cloth when needed.
- Do not clean the unit with substances, such as gasoline, paint thinner or other organic solvents to avoid any damage to the meter.
- Do not clean the LCD or sensor area with water. Lightly wipe with a soft, clean, dry cloth.
- The strip channel must be kept clean. Lightly wipe with a soft, clean, dry cloth each day. Use a damp cloth as needed. Refer to the Maintenance section.
- Follow all local regulations when discarding the unit or its accessories.
- Do not use the unit or the devices outside of the operating temperature ranges: 15-40°C (59-104°F); ≤ 90% RH.

Display	Causes	Solution
	The sensor area is damaged, dirty, or blocked at turn-on, such as a used test strip left in the meter.	Ensure the sensor area is clean and that there are no objects covering the sensor area. Refer to Maintenance. Restart the meter. Contact your local distributor if the sensor area window is broken
	Test strip was removed during the test.	Repeat the test and ensure the test strip remains in place.
	Specimen was applied to the test strip too soon.	Repeat the test and apply specimen after blood drop symbol appears.
	Batteries are discharged but have enough power to run 20 more tests.	Test results will still be accurate, but replace the batteries as soon as possible.
	Batteries are low and meter will not allow more tests until the batteries are replaced.	Replace the batteries then repeat the test.
	Insufficient specimen.	Repeat the test. Apply enough specimen. Use around 25 μL of specimen.
	Expired test strips or incorrect date entered .	Ensure the test strips are within the expiration date printed on the package label. If the test strips are still within the expiration date, check to see if the date was entered correctly. Remove the code chip before starting the setup mode.

	<p>Code chip was removed during testing.</p>	<p>Insert proper code chip. Confirm the code chip matches the test strip code and repeat the test.</p>
	<p>The test strip type does not match the code chip.</p>	<p>Use the proper strip which its type matches the code chip.</p>
	<p>The environment temperature is higher than 40 °C (104°F).</p>	<p>Get the meter in a proper environment where the temperature is between 15-40°C (59-104°F) .</p>
	<p>The environment temperature is lower than 15 °C (59°F)</p>	
<p>CODE</p> 	<p>No code chip in the meter. Code chip is damaged or inserted incorrectly.</p>	<p>Insert the code chip that accompanied the package of test strips. If the code chip is damaged, use a new code chip with the correct code number. If the code chip is inserted incorrectly, remove the code chip and insert it into the code chip slot.</p>

Feature	Specifications
Methodology	Reflectance Photometer
Test Time	≤ 2 min
Measurement Range	TC: 100-500 mg/dL (2.59-12.93 mmol/L, 1 mmol/L=38.66 mg/dL) HDL: 15-120 mg/dL (0.39-3.10 mmol/L, 1 mmol/L=38.66 mg/dL) TG: 45-650 mg/dL (0.51-7.34 mmol/L, 1 mmol/L=88.6 mg/dL)
Specimen	Whole blood, plasma, and serum
Specimen Volume	25 µL
Power Source	Batteries: 3 AAA (LR03)
	Mini USB, 5V dc, 50 mA, use only AC/DC adapter or PC certified and recommended.
Battery Life	85 hours or 1,000 tests
Units of Measurement	mg/dL, mmol/L
Memory	200 records
Automatic Shut Off	5 minutes after last use
Meter Size	142 mm × 70 mm × 27 mm
Display Size	50 mm × 50 mm (1.97" × 1.97")
Weight	130 g (without batteries)
Meter Storage Conditions	0-50°C (32-122°F); ≤ 90% RH; altitude ≤ 2000 m
Operating Conditions	15-40°C (59-104°F); ≤ 90% RH; altitude ≤ 2000 m
Control Strip Storage Conditions	2-30°C (36-86°F); ≤ 90% RH

	Consult instructions for use
REF	Article number
	Manufacturer
LOT	Batch code
	Contains sufficient for <n> tests
STERILE R	Sterilized using irradiation
	Date of manufacture
CE	CE marked
	Do not discard along with household waste
	Keep away from sunlight
	Do not re-use
IVD	In vitro diagnostic medical device
SN	Serial Number
EC REP	Authorized Representative
	Use-by date
	Temperature limit
CODE	Code Number
	Keep Dry
#	Model number
	Importeur

Note: This warranty applies only to the analyzer in the original purchase. It does not apply to the other materials included with the analyzer.

VivaChek Biotech (Hangzhou) Co., Ltd. Warrants to the original purchaser that this meter will be free from defects in materials and workmanship for a period of two years (24 months). The two years starts from the later of the date of original purchase or installation (except as noted below). During the stated two years period, VivaChek shall replace the unit under warranty with a reconditioned unit or, at its option, repair at no charge a unit that is found to be defective. VivaChek shall not be responsible for shipping charges incurred in the repair of such a meter.

This Warranty is subject to the following exception and limitations:

This warranty is limited to repair or replacement due to defects in parts or workmanship. Parts required which were not defective shall be replaced at additional cost. VivaChek shall not be required to make any repairs or replace any parts that are necessitated by abuse, accidents, alteration, misuse, neglect, failure to operate the analyzer in accordance with the operations manual, or maintenance by anyone other than VivaChek. Furthermore, VivaChek assumes no liability from malfunction or damage to meters caused by the use of strips other than strips manufactured by VivaChek. VivaChek reserves the right to make changes in the design of this meter without obligation to incorporate such changes into previously manufactured meters.

Disclaimer of Warranties

This warranty is expressly made in lieu of any and all other warranties express or implied (either in fact or by operation of law) including the warranties of merchantability and fitness for use, which are expressly excluded, and is the only warranty given by VivaChek.

Limitations of Liability

In no event shall VivaChek be liable for indirect, special or consequential damages, even if VivaChek has been advised of the possibility of such damages.

For warranty service, please contact your local distributor.

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Lotus NL B.V.
Hangzhou Bosure Biotech Co., Ltd